


KEY FACTS

Biomarkers for pre-eclampsia management > week 20



The role of angiogenic factors in pre-eclampsia

Although the cause of pre-eclampsia remains unclear, an imbalance of the angiogenic factors Placental Growth Factor (PlGF) and soluble FMS-like Tyrosine Kinase (sFlt-1) seems to induce endothelium dysfunction and therefore may initiate the disease.



Measuring PlGF and sFlt-1 in pregnancy

As separation of curves of healthy vs pre-eclamptic pregnancy is already showing from week 8, PlGF is a perfect biomarker for first trimester pre-eclampsia screening as well as pre-eclampsia diagnosis >20 weeks. Differences in sFlt-1 levels are starting to show after week 20, which makes sFlt-1 the perfect partner for PlGF in the diagnosis of pre-eclampsia > week 20.



● Control ● Women who had pre-eclampsia >5 weeks later ● Women who later had pre-eclampsia ● Women with clinical pre-eclampsia

PlGF belongs to the Vascular Endothelial Growth Factors (VEGF) family, promoting proliferation and survival of endothelial cells.¹

sFlt-1 is a truncated form of the VEGF receptor Flt-1 and binds VEGF and PlGF with high affinity, therefore neutralizing their effects.²

Thermo Scientific™ B·R·A·H·M·S™ PlGF plus KRYPTOR™ is an automated immunofluorescent assay for the quantitative determination of PlGF-1 in human serum and EDTA plasma.

- FAS: 6.7 pg/mL
- Onboard stability: 29 days
- Short incubation time: 29 minutes

Thermo Scientific™ B·R·A·H·M·S™ sFlt-1 KRYPTOR™ is an automated immunofluorescent assay for the quantitative determination of sFlt-1 in human serum and EDTA plasma.


- FAS: 29 pg/mL
- Onboard stability: 29 days
- Short incubation time: 9 minutes

Improved pre-eclampsia diagnosis and short-term prognosis

10% of pregnant women show signs and symptoms of pre-eclampsia. Only **one fifth of them** are actually developing pre-eclampsia.³ Assessment of blood pressure and proteinuria offers only a **poor sensitivity and specificity** in predicting pre-eclampsia.^{4,5} This diagnostic standard is significantly improved by sFit-1 and PIGF serum measurement.


Improved pre-eclampsia diagnosis with sFit-1/PIGF ratio

Women with hypertension due to pre-eclampsia have a significantly higher sFit-1/PIGF ratio than women with any other form of hypertensive disorder.



Short-term prognosis with sFit-1/PIGF ratio in pre-eclamptic women

Women with any subsequent adverse outcome due to pre-eclampsia have a significantly higher sFit-1/PIGF ratio than women with no complications.^{6,9}



Benefits of the sFit-1/PIGF ratio in diagnosing pre-eclampsia

- Significantly higher sensitivity and specificity in diagnosing pre-eclampsia than measuring blood pressure and proteinuria⁷
- Clearly differentiating between pre-eclampsia and other forms of hypertensive disorders⁸

Benefits of the sFit-1/PIGF ratio as a predictor for adverse outcome

- Improved preparation time for delivery if needed^{6,9}
- Improved outbalancing of maternal risk vs fetal prematurity^{6,9}

References 1. Levine et al. N Engl J Med 2004; 350: 672-83 2. Hagmann et al. Clin Chem 2012; 58(5): 837-45 3. Milne et al. BMJ 2009; 339: b3129 4. Verlohren et al. Clin Sci 2012; 122(2): 43-52 5. Zhang et al. Obstet Gyneol 2001; 97: 261-7 6. Rana et al. Circulation 2012; 125 (7): 911-9 7. Hagmann et al. Clin Chem 2012; 58(5): 837-45 8. Verlohren et al. Am J Obstet Gynecol 2010; 202: 161.e1-11 9. Salahuddin et al. Hypert Preg 2016; 35(3): 330-45

Clinical Diagnostics

Thermo Fisher Scientific
B·R·A·H·M·S GmbH
Neuendorfstr. 25
16761 Hennigsdorf
Germany

+49 (0)3302 883 0
+49 (0)3302 883 100 fax
info.brahms@thermofisher.com

Find out more at thermoscientific.com/brahms

© 2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. KRYPTOR is a trademark of CIS bio international, licensed for use by B·R·A·H·M·S, a part of Thermo Fisher Scientific.

Thermo Fisher Scientific products are distributed worldwide; not all intended uses and applications mentioned in this printing are registered in every country.

107431.2